

## Clinical uses of Accent Prime and its status in Korea

You&Chung Skin Clinic, Seoul  
Woonkyong CHUNG

It's been over five and a half years since Accent Prime was introduced in Korea. Accent Prime is a platform device that can be used in a variety of ways by being equipped with various handpieces. Interestingly, while overseas, the ultrasound module is the main attraction, in Korea, the radiofrequency(RF) is more popular.

The most famous RF handpiece of Accent Prime in Korea is the TuneFace. TuneFace is a vacuum-assisted, unipolar RF treatment that is virtually painless and requires no consumables. In addition, the end-point of the treatment can be set differently for each user, so the treatment protocol, i.e. the intensity, time, and number of repetitions, can be customized according to the concept of each clinic.

There are many different types of RF devices that use heat to produce skin tightening and contouring effects. Each device differs in the way it delivers heat, the temperature of the heat, the depth at which the heat is delivered, the amount of heat delivered, and the methods used to reduce pain during the procedure.

So Today, we're going to take a look at what has made TuneFace stand out in the Korean market over the past five and a half years amidst the plethora of RF devices. I'll also share some of the things I've been working on to make TuneFace treatments more effective.